

Figure 1A. Peptide-linker- β chain construct




Figure 1B. Schematic view of peptide-linked MHC binding groove






Figure 1C. Schematic view of soluble peptide-linked MHC-IgG C-region fusion protein




2/58

Figure 2: I-A^d α chain cloning scheme


3/54

Figure 3: I-A^d β chain cloning scheme




4/
56

Figure 3: cont.




T09020-6200000

5/54


Figure 4: HLA-DR1 α chain cloning scheme

6/58

Figure 4: cont.




7/58

Figure 5: HLA-DR1 β chain cloning scheme

8/56


Figure 5: cont.




FOUO 62200600

9/58

Figure 5: cont.




10/58

Figure 6: I-A^S α chain cloning scheme

11/58

Figure 7: I-A^S β chain cloning scheme



12/58

Figure 8 - Oligonucleotides used in constructing MHC vectors

I-A^d/I-A^s PCR primers and cloning oligonucleotides (restriction site are underlined).

OPR100

5'-GGG GGG GCC ATG GCC GAA GAC GAC ATT GAG GCC GAC-3'

OPR101

5'-GGG GCG ACT AGT CCA GTG TTT CAG AAC CGG CTC-3'

OPR107

5'-CCC CCC GAT ATC TCA GGT TCC AGC AGT GGA GAC GAC ATT GAG
GCC G-3'

OPR108

5'-CCC CCC CGG CCG CTA CTT ACG TTT CCA GTG TTT CAG AAC CGG
C-3'

OPR102

5'-GGG GGG GCC ATG GCC GGA AAC TCC GAA AGG CAT TTC G-3'

OPR104

5'-GGG GCG ACT AGT CCA CTC CAC AGT GAT GGG GC-3'

OPR106

5'-CCC CCC CGG CCG TAC CTG AGG ACC ACT CCA CAG TGA TGG-3'

OPR112

5'-CCC CCC GAT ATC ACA GGT GTC TTA AGT GCT AGC GGA GGG GGC
GGA AGC GGC GGA GGG GGA AAC TCC GAA AGG CAT TTC-3'

OPR119

5'-AGC TTG ATA TCA CAG GTG TCT TAA GTG GAG-3'

OPR120-2

5'-CTA GCT CCA CTT AAG ACA CCT GTG ATA TCA-3'

VW310

5'-TCC GGA GGC GGC GGA GAC TCC GAA AGG CAT TTC G-3'

VW309

5'-CGA TCG CTA GCG GCG GTG GTG GTT CCG GTG GCG GCG GAG-3'

OPR136

5'-CCC CCC AGG CTT CCC GGG CCA CCA TGC CTT GCA GCA GAG CTC
TG-3'

OPR139

5'-CCC CCC GAG CTC GAA TTC TCA TAA AGG CCC TGG GTG TCT G-3'

FOUO 20060229 10504

13/59

Figure 8 - cont.

OPR132

5'-CCC CCC AAG CTT CCC GGG CCA CCA TGG CTC TGC AGA TCC CCA
GC-3'

OPR133

5'-CCC CCC ACT TAA GGT CCT TGG GCT GCT CAG CAC C-3'

OPR134

5'-CCC CCC CCA TCA CTG TGG AGT GGA GGG-3'

OPR135

5'-CCC CCC GAG CTC GAA TTC TCA CTG CAG GAG CCC TGC TGG-3'

HLA-DR1 PCR primers and cloning oligonucleotides.

DR1A-F

5'-GGG GGG AAG CTT ATG ATC AAA GAA GAA CAT GTG ATC ATC-3'

DR1A-B

5'-GCG GCG GGA TCC GTT CTC TGT AGT CTC TGG GAG AGG-3'

DR1B-F

5'-GGG GGG AAG CTT ATG GGG GAC ACC CCA CCA CGT TTC TTG TGG
CAG C-3'

AF-N

5'-GGG GGG GCC ATG GCC ATC AAA GAA GAA CAT GTG ATC ATC-3'

AB-S

5'-GCG GCG ACT AGT GTT CTC TGT AGT CTC TGG GAG AGG-3'

OPR124

5'-GGG GGG AAG CTT GAT ATC TCA GCT TCC AGC AGT AGT ATC AAA
GAA GAA CAT GTG ATC-3'

OPR125

5'-GGG GGG CGG CCG CTA CTT ACC TTT CTC TGG GAG AGG GCT TGG
AGC-3'

DR1B-B

5'-GCG GCG GGA TCC CTT GCT CTG TGC AGA TTC AGA CC-3'

BF-NN

5'-GGG GGG GCC ATG GCC GCA TCC GCT ACC GGG GAC ACC CCA CCA
CGT TTC TTG-3'

F09020-0200000

14/58

BB-5

088121

OPR:22

CPA:23

Peptide oligonucleotides.

Ova 323-332

OPR:10

000000

Ova EGIR

Quesada

OPR.116

Ova A532Y

OPR117

OPR115

5'-CTA GCA CGA CCA GGT TCG TTG ATT TCA TAG TGA GCA GCG TGA
ACA GCG TGA GAG ATA C-3'

15/58

Figure 3 - cont.

HEL 74-86

OPR140

5'-TTA AGT AAC CTG TGC AAC ATC CCC TGC AGC GCC CTG CTG AGC
TCC G-3'

OPR141

5'-CTA GCG GAG CTC AGC AGG GCG CTG CAGBGGG ATG TTG CAC AGG
TTA C-3'

NP 404-415

OPR128

5'-TTA AGT CAG ATC AGC GTG CAG CCC GCC TTC AGC GTG CAG G-3'

OPR129

5'-CTA GCG TGC ACG CTG AAG GCG GGC TGA ACG CTG ATC TGA C-3'

HA 307-319

OPR130

5'-TTA AGT CCC AAG TAC GTG AAG CAG AAC ACC CTG AAG CTG GCC
ACC G-3'

OPR131

5'-CTA GCG GTG GCC AGC TTC AGG GTG TTC TGC TTC ACG TAC TTG
GGA C-3'

MBP 91-103

VW315

5'-TTA AGT CAC TAT GGC TCC CTG CCG CAG AAG TCC CAG CAC GGG
CGC G-3'

VW316

5'-CTA GCG CGC CCG TGC TGG GAC TTC TGC GGC AGG GAG CCA TAG
TGA C-3'

PLP 139-151

VW313

5'-TTA CAG CAC TCC CTG GGC AAG TGG CTG GGC CAC CCG GAC AAG
TTC G-3'

VW314

5'-CTA GCG AAC TTG TTC GGG TGG CCC AGC CAC TTG CCC AGG GAG
TGA C-3'

FO96020-6200000

Figure 8 - cont.

16/58

MBP 1-14

VW317

5' - TTA AGT ATG GCA TCC CAG AAG CGC CCG TCC CAG CGC TCC AAG
TAC CTG G-3'

VW316

5' - CTA GCC AGG TAC TTG GAG CGC TGG GAC GGG CGC TTC TGG GAT
GCC ATA C-3'

T09020 62200000


Figure 9A. Soluble I-A^d α chain construct

Figure 9B.

Soluble I-A^d β chain construct

Restriction sites for insertion of
oligonucleotides encoding peptides of interest




Figure 9C.


Soluble I-A^s α chain construct

Figure 9 - cont.
 Figure 9D. Soluble I-A^s β chain construct




Figure 9E. Soluble HLA-DR1 α chain construct







Figure 9F. Soluble HLA-DR1 β chain construct



19/58

Figure 10. Original Mammalian Cell Expression Vectors



20/58

Figure 11. The 27 KB Kappa and the 17 Kb Heavy Chain EcoR V and Eag I Mutated Constructs

Figure 11A






Figure 11B

[illegible]


21/59

Figure 12. PCR Site Directed Mutagenesis for Introducing EcoRV and EagI Restriction Sites into Kappa Chain 2.7kb Insert



22/58

Figure 13. STRATEGY FOR CONSTRUCTING MEC CLASS II α / KAPPA CONSTANT GENE IN MAMMALIAN CELL EXPRESSION VECTOR




23/58

Figure 14. Primers used for Sequencing Mutated 2.7 Kb Fragment

<u>Primer List</u>	<u>Sequence</u>
PMC-33	[5GCTCAGCTGTCTTGTTTCAGTACTGATG3']
PMC-77	[5GTAAGTAGCGGCCG3']
PMC-111	[5GGTATGTAAAAATAAACATCACAG3']
PMC-114	[5GCTTTGCTTACGGAGTTACTG3']

FOUO 62200000

24/58

Figure 15. Strategy for M13 Mutagenesis and Cloning of the MEC II β Variable Gene

$$25/58$$

Figure 16A






Figure 16B

[illegible]


26/58

Figure 17 - Full length peptide linked MHC expression vectors



27/58

Figure 17 - cont.



28/58

Figure 17 - cont.






Figure 18A. Full length I-Ad α chain insertFigure 18B. Full length I-Ad β chain insert

FIGURE 19 (Sheet 1 of 7)

Full-length peptide linked MHC expression vectors





FIGURE 19 (Sheet 3 of 7)



33/58

FIGURE 19 (Sheet 4 of 7)






FIGURE 19 (Sheet 4 of 7)

35/59
FIGURE 19 (Sheet 6 of 7)

36/58
FIGURE 19 (Sheet 7 of 7)



37/58
FIGURE 20

Oligonucleotides used in cloning

OPR132

I-A^d β signal peptide front primer with Kozak consensus for CellTech vector -
HindIII/XmaI sites

5'-CCC CCC AAG CTT CCC GGG CCA CCA TGG CTC TGC AGA TCC CCA
GC-3'

OPR133

I-A^d β signal peptide back primer with Kozak consensus for CellTech vector - AflII site

5'-CCC CCC ACT TAA GGT CCT TGG GCT GCT CAG CAC C-3'

OPR134

I-A^d β transmembrane front primer for CellTech vector - BstXI sites

5'-CCC CCC CCA TCA CTG TGG AGT GGA GGG-3'

OPR135

I-A^d β transmembrane back primer for CellTech vector - SstI, EcoRI sites

5'-CCC CCC GAG CTC GAA TCC TCA CTG CAG GAG CCC TGC TGG-3'

OPR136

I-A^d α signal peptide front primer with Kozak consensus for CellTech vector -
HindIII/XmaI sites

5'-CCC CCC AAG CTT CCC GGG CCA CCA TGC CGT GCA GCA GAG CTC
TG-3'

OPR139

I-A^d α transmembrane primer for CellTech vector - SstI/EcoRI sites

5'-CCC CCC GAG CTC GAA TCC TCA TAA AGG CCC TGG GTG TCT G-3'

B7-1-2F

Murine B7-1 front primer with Kozak consensus for CloneTech vector - NotI site

5'-CCC CCC CCG CGG CCG CCC CAC CAT GGG ACT GAG TAA CAT TCT
C-3'

B7-1-2B

Murine B7-1 BACK primer for CloneTech vector - NotI site


5'-CCC CCC GCG GCC GCT TTA AAA ACA TGT ATC ACT TTT-3'

T09040-620000

33/58


FIGURE 21


NSO/ClassII/OVA clones Stimulate IL-2 Production from DO11.10



39/58
FIGURE 22


INTRAMUSCULAR IAd/OVA & B7 DNA INJECTIONS PROLIFERATION ASSAY




40/58
FIGURE 23**INTRADERMAL IAd/OVA & IAd/HEL DNA INJECTIONS
PROLIFERATION ASSAY
4, 7 & 14 DAYS POST INJECTION**


41/58


FIGURE 24



090029-0200


42/56
FIGURE 25 (SHEET 1 OF 4)

43/58
FIGURE 25 (SHEET 2 OF 4)

44/54
FIGURE 25 (SHEET 3 OF 4)

45/54

FIGURE 25 (SHEET 4 OF 4)



46/58
FIGURE 26

JLA-005

5'-CCCCCCCCTATGGCCGCTAGCGGAGGGGGCGGAAGC-3'

JLA-007

5'-CCCGGGGCTCGAGTGAAGACGACATTGAGGCCGAC-3'

JLA-009

5'-CCCCCACTAGTCCACTCCACAGTGATGGGGCT-3'

JLA-010

5'-CCCCCCCCCGGGACCAGTGTTTCAGAACCGGCTCCTC-3'

JLA-301

5'-TCGAGGAACCGCCACCGCCAGAACCGCCGCCACCGGA-
ACCACCACCGCCGCTGCCACCGCCACCA-3'

JLA-302

5'-CTAGTGGTGGCGGTGGCAGCGGCGGTGGTGGTTCCGG-
TGGCGGCGGTCTGGCGGTGGCGGTTC-3'

OPR-142

5'-CTTGGAATCTTGACTAAGAGG-3'

JS-305

5'-CAGGTGGAATTCTCATTCCATCGGCATGTACTCTTCTT-
CCTCCCAGTGTTTCAGAACCGG-3'

T09020 020000

47/54
FIGURE 27 (SHEET 1 OF 4)

```

      10      20      30      40      50
      .      .      .      .      .
CCACC ATG GCT CTG CAG ATC CCC AGC CTC CTC CTC TCA GCT GCT GTG GTG
GGTGG TAC CGA GAC GTC TAG GGG TCG GAG GAG GAG AGT CGA CGA CAC CAC
      M   A   L   Q   I   P   S   L   L   L   S   A   A   V   V>
<----- I-Ad β chain leader ----->

      60      70      80      90
      .      .      .      .
GTG CTG ATG GTG CTG AGC AGC CCA AGG ACC TTA AGT ATC TCT CAG GCT
CAC GAC TAC CAC GAC TCG TCG GGT TCG AAT TCA TAG AGA GTC CGA
      V   L   M   V   L   S   S   P   R   T   L   S   I   S   Q   A>
----->-----<-----

    100      110      120      130      140
    .      .      .      .      .
GTT CAC GCT GCT CAC GCT GAA ATC AAC GAA GCT GGT CGT GCT AGC GGA
CAA GTG CGA CGA GTG CGA CTT TAG TTG CTT CGA CCA GCA CGA TCG CCT
      V   H   A   A   H   A   E   I   N   E   A   G   R   A   S   G>
----- OVA 323-339 ----->-----<-----

    150      160      170      180      190
    .      .      .      .      .
GGG GGC GGA AGC GGC GGA GGG CGA AAC TCC GAA AGG CAT TTC GTG GTC
CCC CCG CCT TCG CCG CCT CCC CCT TTG AGG CTT TCC GTA AAG CAC CAG
      G   G   G   S   G   G   G   G   N   S   E   R   H   F   V   V>
-- 10 amino acid linker -->----- I-Ad β-1 domain -----<-----

    200      210      220      230      240
    .      .      .      .      .
CAG TTC AAG GGC GAG TGC TAC TAC ACC AAC GGG ACG CAG CGC ATA CGG
GTC AAG TTC CCG CTC ACC ATG ATG TCG TTG CCC TGC GTC GCG TAT GCC
      Q   F   K   G   E   C   Y   Y   T   N   G   T   Q   R   I   R>
----->-----<-----

    250      260      270      280      290
    .      .      .      .      .
CTC GTG ACC AGA TAC ATC TAC AAC CGG GAG GAG TAC GTG CGC TAC GAC
GAG CAC TCG TCT ATG TAG ATG TTG GCC CTC CTC ATG CAC GCG ATG CTG
      L   V   T   R   Y   I   Y   N   R   E   E   Y   V   R   Y   D>
----->-----<-----

    300      310      320      330
    .      .      .      .
AGC GAC GTG GGC GAG TAC CGC GCG GTG ACC GAG CTG GGG CGG CCA GAC
TCG CTC CAC CCG CTC ATG GCG GCG CAC TCG CTC GAC CCC GCC GGT CTG
      S   D   V   G   E   Y   R   A   V   T   E   L   G   R   P   D>
----->-----<-----

    340      350      360      370      380
    .      .      .      .      .
GCC GAG TAC TCG AAC AGC CAG CCG GAG ATC CTG GAG CGA ACG CGG GCC
CGG CTC ATG ACC TTG TCG GTC GGC CTC TAG GAC CTC GCT TGC GCC CGG
      A   E   Y   W   N   S   Q   P   E   I   L   E   R   T   R   A>
----->-----<-----

```

T03020-020000

48/58
FIGURE 27 (SHEET 2 OF 4)

```

390      400      410      420      430
GAG GTG GAC ACG GCG TGC AGA CAC AAC TAC GAG GGG CCG GAG ACC AGC
CTC CAC CTG TGC CGC ACG TCT GTG TTG ATG CTC CCC GGC CTC TGG TCG
E V D T A C R H N Y E G P E T S>
-----

440      450      460      470      480
ACC TCC CTG CCG CGG CTT GAA CAG CCC AAT GTC GCC ATC TCC CTG TCC
TGG AGG GAC GCC GCC GAA CTT GTC GGG TTA CAG CGG TAG AGG GAC AGG
T S L R R L E Q P N V A I S L S>
--- I-Ad β-1 domain ->----- I-Ad β-2 domain -----

490      500      510      520      530
AGG ACA GAG GCC CTC AAC CAC CAC AAC ACT CTG GTC TGT TCG GTG ACA
TCC TGT CTC CGG GAG TTG GTG GTC TTG TGA GAC CAG ACA AGC CAC TGT
R T E A L N H H N T L V C S V T>
-----

540      550      560      570
GAT TTC TAC CCA GCC AAG ATC AAA GTG CGC TGG TTC AGG AAT GGC CAG
CTA AAG ATG GGT CGG TTC TAG TTT CAC GCG ACC AAG TCC TTA CCG GTC
D F Y P A K I K V R W F R N G Q>
-----

580      590      600      610      620
GAG GAG ACA GTG GGG GTC TCA TCC ACA CAG CTT ATT AGG AAT GGG GAC
CTC CTC TGT CAC CCC CAG AGT AGG TGT GTC GAA TAA TCC TTA CCC CTG
E E T V G V S S T Q L I R N G D>
-----

630      640      650      660      670
TGG ACC TTC CAG GTC CTG GTC ATG CTG GAC ATG ACC CCT CAT CAG GGA
ACC TGG AAG GTC CAG GAC CAG TAC GAC CTC TAC TGG GGA GTA GTC CCT
W T F Q V L V M L E M T P H Q G>
-----

680      690      700      710      720
GAG GTC TAC ACC TGC CAT GTG GAG CAT CCC AGC CTC AAG AGC CCC ATC
CTC CAG ATG TGG ACG GTA CAC CTC GTA GGG TCG GAC TTC TCG GGG TAG
E V Y T C H V E H P S L R S P I>
----- I-Ad β-2 domain -----

730      740      750      760      770
ACT GTG GAG TGG ACT AGT GGT GGC GGT AGC GGC GGT GGT GGT TCC
TGA CAC CTC ACC TGA TCA CCA CCG CCA CCG TCG CCG CCA CCA CCA AGG
T V E W T S G G G G S G G G G S>
-----><----- 24 amino acid linker -----

```

T03020-020000

49/59
FIGURE 27 (SHEET 3 OF 4)

```

      780      790      800      810
      .      .      .      .
GGT GGC GGC GGT TCT GGC GGT GGC GGT TCC TCC AGT GAA GAC GAC ATT
CCA CCG CCG CCA AGA CCG CCA AGG AGG TCA CTT CTG CTG TAA
G G G G S G G G G S S S E D D I>
----->-----

820      830      840      850      860
      .      .      .      .      .
GAG GCC GAC CAC GTA GGC TTC TAT GGT ACA ACT GTT TAT CAG TCT CCT
CTC CCG CTG GTG CAT CCG AAG ATA CCA TGT TGA CAA ATA GTC AGA GGA
E A D H V G F Y G T T V Y Q S P>
----- I-Ad α-1 domain -----

870      880      890      900      910
      .      .      .      .      .
GGA GAC ATT GCC CAG TAC ACA CAT GAA TTT GAT GGT GAT GAG TTG TTC
CCT CTG TAA CCG CTC ATG TGT GTA CTT AAA CTA CCA CTA CTC AAC AAG
G D I G Q Y T H E F D G D E L F>
-----

920      930      940      950      960
      .      .      .      .      .
TAT GTG GAC TTG GAT AAG AAG AAA ACT GTC TGG AGG CTT CCT GAG TTT
ATA CAC CTG AAC CTA TTC TTC TTT TGA CAG ACC TCC GAA GGA CTC AAA
Y V D L D K K K T V W R L F E F>
-----

970      980      990      1000      1010
      .      .      .      .      .
GGC CAA TTG ATA CTC TTT GAG CCC CAA GGT GGA CTG CAA AAC ATA GCT
CCG GTT AAC TAT GAG AAA CTC GGG GTT CCA CCT GAC GTT TTG TAT CGA
G Q L I L F E P Q G G L Q N I A>
-----

1020      1030      1040      1050
      .      .      .      .
GCA GAA AAA CAC AAC TTG GGA ATC TTG ACT AAG AGG TCA AAT TTC ACC
CGT CTT TTT GTG TTG AAC CCT TAG AAC TGA TTC TCC AGT TTA AAG TGG
A E K H N L G I L T K R S N F T>
----- I-Ad α-1 domain -----

1060      1070      1080      1090      1100
      .      .      .      .      .
CCA GCT ACC AAT GAG GCT CCT CAA GCG ACT GTG TTC CCC AAG TCC CCT
GGT CGA TGG TTA CTC CGA CGA GTT CGC TGA CAC AAG GGG TTC AGG GGA
P A T N E A P Q A T V F P K S P>
----->----- I-Ad α-2 domain -----

1110      1120      1130      1140      1150
      .      .      .      .      .
GTG CTG CTG GGT CAG CCC AAC ACC CTT ATC TGC TTT GTG GAC AAC ATC
CAC GAC GAC CCA GTC GGG TTG TGG GAA TAG ACG AAA CAC CTG TTG TAG
V L L G Q P N T L I C F V D N I>
-----

```

T09020-62E00000

5/54
FIGURE 28 (SHEET 1 OF 4)

```

      10           20           30           40           50
      .           .           .           .           .
CCACC ATG GCT CTG CAG ATC CCC AGC CTC CTC CTC TCA GCT GCT GTG GTG
GGTGG TAC CGA GAC GTC TAG GGG TCG GAG GAG GAG AGT CGA CGA CAC CAC
      M   A   L   Q   I   P   S   L   L   L   S   A   A   V   V>
<----- I-Ad β chain leader ----->

      60           70           80           90
      .           .           .           .
GTC CTG ATG GTG CTG AGC AGC CCA AGG ACC TTA AGT ATC TCT CAG GCT
CAC GAC TAC CAC GAC TCG TCG GGT TCC TCG AAT TCA TAG AGA GTC CGA
      V   L   M   V   L   S   S   P   R   T   L   S   I   S   Q   A>
----->-----<----->-----

    100           110           120           130           140
    .           .           .           .           .
GTT CAC GCT GCT CAC GCT GAA ATC AAC GAA GCT GGT CGT GCT AGC GGA
CAA GTG CGA CGA GTG CGA CTT TAG TTG CTT CGA CCA GCA CGA TCG CCT
      V   H   A   A   H   A   E   I   N   E   A   G   R   A   S   G>
----- OVA 323-339 ----->-----<----->-----

    150           160           170           180           190
    .           .           .           .           .
GGG GGC GGA AGC GGC GGA GGG GGA AAC TCG GAA AGG CAT TTC GTG GTC
CCC CCG CCT TCG CCG CCT CCC CCT TTG AGG CTT TCC GTA AAG CAC CAG
      G   G   G   S   G   G   G   G   N   S   E   R   H   F   V   V>
-- 10 amino acid linker --><----- I-Ad β-1 domain ----->-----

    200           210           220           230           240
    .           .           .           .           .
CAG TTC AAG GGC GAG TGC TAC TAC ACC AAC GGG ACG CAG CGC ATA CGG
GTC AAG TTC CCG CTC ACG ATG ATG TGG TTG CCC TCC GTC GCG TAT GCC
      Q   F   K   G   E   C   Y   Y   T   N   G   T   Q   R   I   R>
----->-----<----->-----

    250           260           270           280           290
    .           .           .           .           .
CTC GTG ACC AGA TAC ATC TAC AAC CGG GAG GAG TAC GTG CGC TAC GAC
GAG CAC TGG TCT ATG TAG ATG TTG GCC CTC CTC ATG CAC GCG ATG CTG
      L   V   T   R   Y   I   Y   N   R   E   E   Y   V   R   Y   D>
----->-----<----->-----

    300           310           320           330
    .           .           .           .
AGC GAC GTG GGC GAG TAC CSC GCG GTG ACC GAG CTG GGG CGG CCA GAC
TCG CTG CAC CCG CTC ATG GCG CGC CAC TGG CTC GAC CCC GCC GGT CTG
      S   D   V   G   E   Y   R   A   V   T   E   L   G   R   P   D>
----->-----<----->-----

    340           350           360           370           380
    .           .           .           .           .
GCC GAG TAC TGG AAC AGC CAG CCG GAG ATC CTG GAG CGA ACG CGG GCC
CGG CTC ATG ACC TTG TCG GTC GCG CTC TAG GAC CTC GCT TGC GCC CGG
      A   E   Y   W   N   S   Q   P   E   I   L   E   R   T   R   A>
----->-----<----->-----

```

T03020-020000

52/58
FIGURE 28 (SHEET 2 OF 4)

```

390      400      410      420      430
GAG GTG GAC ACG GCG TGC AGA CAC AAC TAC GAG GGG CCG GAG ACC AGC
CTC CAC CTG TGC CGC ACG TCT GTG TTG ATG CTC CCC GGC CTC TGG TCG
E V D T A C R H N Y E G P E T S>
-----

440      450      460      470      480
ACC TCC CTG CCG CGG CTT GAA CAG CCC AAT GTC GCC ATC TCC CTG TCC
TGG AGG GAC GCC GGC GAA CTT GTC GGG TTA CAG CCG TAG AGG GAC AGG
T S L R R L E Q P N V A I S L S>
--- I-Ad β-1 domain -->----- I-Ad β-2 domain -----

490      500      510      520      530
AGG ACA GAG GCC CTC AAC CAC CAC AAC ACT CTG GTC TGT TCG GTG ACA
TCC TGT CTC CGG GAG TTG GTG GTG TTG TGA GAC CAG ACA AGC CAC TGT
R T E A L N H H N T L V C S V T>
-----

540      550      560      570
GAT TTC TAC CCA GCC AAG ATC AAA GTG CGC TGG TTC AGG AAT GGC CAG
CTA AAG ATG GGT CGG TTC TAG TTT CAC GCG ACC AAG TCC TTA CCG GTC
D F Y P A K I K V R W F R N G Q>
-----

580      590      600      610      620
GAG GAG ACA GTG GGG GTC TCA TCC ACA CAG CTT ATT AGG AAT GGG GAC
CTC CTC TGT CAC CCC CAG ACT AGG TGT GTC GAA TAA TCC TTA CCC CTG
E E T V G V S S T Q L I R N G D>
-----

630      640      650      660      670
TGG ACC TTC CAG GTC CTG GTC ATG CTG GAG ATG ACC CCT CAT CAG GGA
ACC TGG AAG GTC CAG GAC CAG TAC GAC CTC TAC TGG GGA GTA GTC CCT
W T F Q V L V M L E M T P H Q G>
-----

680      690      700      710      720
GAG GTC TAC ACC TGC CAT GTG GAG CAT CCC ACC CTG AAG AGC CCC ATC
CTC CAG ATG TGG ACC GTA CAC CTC GTA GGG TCG GAC TTC TCG GGC TAG
E V Y T C H V E H P S L K S P I>
----- I-Ad β-2 domain -----

730      740      750      760      770
ACT GTG CAG TGG ACT ACT GGT GCG GGT GCG AGC GCG GGT GGT GGT TCC
TGA CAC CTC ACC TGA TCA CCG CCA CCG TCG CCG CCA CCA CCA AGG
T V E W T S G G G S G G G S>
-----><----- 24 amino acid linker -----

```

F03020 C260000

53/58
FIGURE 28 (SHEET 3 OF 4)

780 790 800 810
 GGT GGC GGC GGT TCT GGC GGT GGC TCC TCG AGT GAA GAC GAC ATT
 CCA CCG CCG CCA AGA CCG CCA CCG CCA AGG AGC TCA CTT CTG CTG TAA
 G G G G S G G G S S S E D D I>
 -----><-----

820 830 840 850 860
 GAG GCC GAC CAC GTA GGC TTC TAT GGT ACA ACT GTT TAT CAG TCT CCT
 CTC CCG CTG GTG CAT CCG AAG ATA CCA TGT TGA CAA ATA GTC AGA GGA
 E A D H V G F Y G T T V Y Q S P>
 ----- I-Ad α -1 domain -----

870 880 890 900 910
 GGA GAC ATT GGC CAG TAC ACA CAT GAA TTT GAT GGT GAT GAG TTG TTC
 CCT CTG TAA CCG GTC ATG TGT GTA CTT AAA CTA CCA CTA CTC AAC AAG
 G D I G Q Y T H E F D G D E L F>

920 930 940 950 960
 TAT GTG GAC TTG GAT AAG AAG AAA ACT GTC TCG AGG CTT CCT GAG TTT
 ATA CAC CTG AAC CTA TTC TTC TTT TGA CAG ACC TCC GAA GGA CTC AAA
 Y V D L D K K K T V W R L P E F>

970 980 990 1000 1010
 GGC CAA TTG ATA CTC TTT GAG CCC CAA GGT GGA CTG CAA AAC ATA GGT
 CCG GTT AAC TAT GAG AAA CTC GGG GTT CCA CCT GAC GTT TTG TAT CGA
 G Q L I L F E P Q G G L Q N I A>

1020 1030 1040 1050
 GCA GAA AAA CAC AAC TTG GGA ATC TTG ACT AAG AGG TCA AAT TTC ACC
 CGT CTT TTT GTG TTG AAC CCT TAG AAC TGA TTC TCC AGT TTA AAG TGG
 A E K H N L G I L T K R S N F T>
 ----- I-Ad α -1 domain -----

1060 1070 1080 1090 1100
 CCA GCT ACC AAT GAG GCT CCT CAA GCG ACT GTG TTC CCC AAG TCC CCT
 GGT CGA TCG TTA CTC CGA GGA GTT CCG TGA CAC AAG GCG TTC AGG GGA
 P A T N E A P Q A T V F P K S P>
 -----><----- I-Ad α -2 domain -----

1110 1120 1130 1140 1150
 GTG CTG CTG GGT CAG CCC AAC ACC CTT ATC TGC TTT GTG GAC AAC ATC
 CAC GAC GAC CCA GTC GCG TTG TCG GAA TAG ACG AAA CAC CTG TTG TAG
 V L L G Q P N T L I C F V D N I>

T03020-020000

54/58

FIGURE 28 (SHEET 4 OF 4)

```

1160      1170      1180      1190      1200
TTC CCA CCT GTG ATC AAC ATC ACA TGG CTC AGA AAT AAG AAG TCA GT
AAG GGT GGA CAC TAG TTC TAG TGT ACC GAG TCT TTA TCG TTC ACT CAG
F P P V I N I T W L R N S K S V>
-----
1210      1220      1230      1240      1250
ACA GAC GGC GTT TAT GAG ACC AGC TTC CTC GTC AAC CGT GAC CAT TCC
TGT CTG CCG CAA ATA CTC TGG TCG AAG GAG CAG TTG GCA CTG GTA AGG
T D G V Y E T S F L V N R D H S>
-----
1260      1270      1280      1290
TTC CAC AAG CTG TCT TAT CTC ACC TTC ATC CCT TCT GAT GAT GAC ATT
AAG GTG TTC GAC AGA ATA GAG TGG AAG TAG GGA AGA CTA CTA CTG TAA
F H K L S Y L T F I P S D D D I>
-----
1300      1310      1320      1330      1340
TAT GAC TGC AAG GTG GAG CAC TGG GGC CTG GAG GAG CCG GTT CTG AAA
ATA CTG ACG TTC CAC CTC GTG ACC CCG GAC CTC CTC GGC CAA GAC TTT
Y D C K V E H W G L E E P V L K>
-----
I-Ad  $\alpha$ -2 domain -----
1350      1360      1370      1380      1390
CAC TGG GAA CCT GAG ATT CCA GGC CCC ATG TCA GAG CTG ACA GAA ACT
GTG ACC CTT GGA CTC TAA GGT CCG GCG TAC AGT CTC GAC TGT CTT TGA
H N E P E I P A P M S E L T E T>
-----
<----- I-Ad  $\alpha$ -TM domain -----
1400      1410      1420      1430      1440
GTG GTG TGT GCC CTG GCG TTG TCT GTG GGC CTT GTG GGC ATC GTG GTG
CAC CAC ACA CCG GAC CCC AAC AGA CAC CCG GAA CAC CCG TAG CAC CAC
V V C A L G L S V G L V G I V V>
-----
1450      1460      1470      1480      1490
GGC ACC ATC TTC ATC ATT CAA GGC CTG CGA TCA GGT GGC ACC TCC AGA
CCG TGG TAG AAG TAG TAA GTT CCG GAC GGT AGT CCA CCG TCG AGG TCT
G T I F I I Q G L R S G G T S R>
-----
1500
CAC CCA GGC CCT TTA TGA
GTG GGT CCC GGA AAT ACT
H P G P L *>
- I-Ad  $\alpha$ -TM domain ->

```

T03020-620000

55/59
FIGURE 29 (SHEET 1 OF 4)

```

      10      20      30      40      50
      .      .      .      .      .
CCACC ATG GCT CTG CAG ATC CCC AGC CTC CTC CTC TCA GCT GCT GTG GTG
GGTGG TAC CGA GAC GTC TAG GGG TCG CAG GAG GAG AGT CGA CGA CAC CAC
      M  A  L  Q  I  P  S  L  L  L  S  A  A  V  V>
<----- I-Ad β chain leader ----->

      60      70      80      90
      .      .      .      .
GTG CTG ATG GTG CTG AGC AGC CCA AGG ACC TTA AGT ATC TCT CAG GCT
CAC GAC TAC CAC GAC TCG TCG GGT TCC TGG AAT TCA TAG AGA GTC CGA
      V  L  M  V  L  S  S  P  R  T  L  S  I  S  Q  A>
----->-----<-----

    100      110      120      130      140
    .      .      .      .      .
GTT CAC GCT GCT CAC GGT GAA ATC AAC GAA GCT GGT CGT GCT AGC GGA
CAA GTG CGA CGA GTG CGA CTT TAG TTG CTT CGA CCA GCA CGA TCG CCT
      V  H  A  A  H  A  E  I  N  E  A  G  R  A  S  G>
----- OVA 323-339 ----->-----<-----

    150      160      170      180      190
    .      .      .      .      .
GGG GGC GGA AGC GGC GGA GGG GGA AAC TCC GAA AGG CAT TTC GTG GTC
CCC CCG CCT TCG CCG CCT CCC CCT TTG AGC CTT TCC GTA AAG CAC CAG
      G  G  G  S  G  G  G  G  N  S  E  R  H  F  V  V>
-- 10 amino acid linker -->----- I-Ad β-1 domain ----->-----<-----

    200      210      220      230      240
    .      .      .      .      .
CAG TTC AAG GGC GAG TGC TAC TAC ACC AAC GGG ACC CAG CGC ATA CGG
GTC AAG TTC CCG CTC ACG ATG ATG TGG TTG CCC TGC GTC GCG TAT GCC
      Q  F  K  G  E  C  Y  Y  T  N  G  T  Q  R  I  R>
----->-----<-----

    250      260      270      280      290
    .      .      .      .      .
CTC GTG ACC AGA TAC ATC TAC AAC CCG GAG GAG TAC GTG CGC TAC GAC
GAG CAC TGG TCT ATG TAG ATG TTG GCC CTC CTC ATG CAC GCG ATG CTG
      L  V  T  R  Y  I  Y  N  R  E  E  Y  V  R  Y  D>
----->-----<-----

    300      310      320      330
    .      .      .      .
AGC GAC GTG GGC GAG TAC CCG CCG GTG ACC GAG CTC GGG CCG CCA GAC
TCG CTC CAC CCG CTC ATG GCG CCG CAC TGG CTC GAC CCC GCC GGT CTC
      S  D  V  G  E  Y  R  A  V  T  E  L  G  R  P  D>
----->-----<-----

    340      350      360      370      380
    .      .      .      .      .
GCC GAG TAC TGG AAC AGC CAG CCG GAG ATC CTC GAG CGA ACG CCG GCC
CGG CTC ATG ACC TTG TCG GTC GCG CTC TAG GAC CTC GCT TGC GCC CCG
      A  E  Y  W  N  S  Q  P  E  I  L  E  R  T  R  A>
----->-----<-----

```

FIGURE 29 (SHEET 1 OF 4)

56/58
FIGURE 29 (SHEET 2 OF 4)

```

390      400      410      420      430
GAG GTG GAC ACG GCG TGC AGA CAC AAC TAC GAG GGG CCG GAG ACC AGC
CTC CAC CTG TGC CGC ACG TCT GTG TTG ATG CTC CCC GGC CTC TGG TCG
E V D T A C R H N Y E G P E T S>
-----

440      450      460      470      480
ACC TCC CTG CCG CGG CTT GAA CAG CCC AAT GTC GCC ATC TCC CTG TCC
TGG AGG GAC GCC GCC GAA CTT GTG GGG TTA CAG CGG TAG AGG GAC AGG
T S L R R L E Q P N V A I S L S>
--- I-Ad β-1 domain ->----- I-Ad β-2 domain -----

490      500      510      520      530
AGG ACA GAG GCC CTC AAC CAC CAC AAC ACT CTG GTC TGT TCG GTG ACA
TCC TGT CTC CGG GAG TTG GTG GTG TTG TGA GAC CAG ACA AGC CAC TGT
R T E A L N H H N T L V C S V T>
-----

540      550      560      570
GAT TTC TAC CCA GCC AAG ATC AAA GTG CGC TGG TTC AGG AAT GGC CAG
CTA AAG ATG GGT CGG TTC TAG TTT CAC GCG ACC AAG TCC TTA CCG GTC
D F Y P A K I K V R W F R N G Q>
-----

580      590      600      610      620
GAG GAG ACA GTG GGG GTC TCA TCC ACA CAG CTT ATT AGG AAT GGG GAC
CTC CTC TGT CAC CCC CAG AGT AGG TGT GTC GAA TAA TCC TTA CCC CTC
E E T V G V S S T Q L I R N G D>
-----

630      640      650      660      670
TGG ACC TTC CAG GTC CTC GTC ATG CTC GAG ATG ACC CCT CAT CAG GGA
ACC TGG AAG GTC CAG GAC CAG TAC GAC CTC TAC TGG GGA GTA GTC CCT
W T F Q V L V M L E M T P H Q G>
-----

680      690      700      710      720
GAG GTC TAC ACC TGC CAT GTG GAG CAT CCC AGC CTG AAG AGC CCC ATC
CTC CAG ATG TGG ACG GTA CAC CTC GTA GGG TCG GAC TTC TCG GCG TAG
E V Y T C H V E H P S L K S P I>
----- I-Ad β-2 domain -----

730      740      750      760      770
ACT GTG GAG TGG ACT AGT GGT GGC GGT GGC AGC GGC GGT GGT GGT TCC
TGA CAC CTC ACC TGA TCA CCA CCG CCA CCG TCG CCG CCA CCA CCA AGG
T V E W T S G G G G S G G G G S>
-----><----- 24 amino acid linker -----

```

T03020-6200000

FIGURE 29 (SHEET 3 OF 4)

```

      780      790      800      810
      .      .      .      .
GCT GGC GGC GGT TCT GGC GGT GGC GGT TCC TCG AGT GAA GAC GAC ATT
CCA CCG CCG CCA AGA CCG CCA CCG CCA AGC AGC TCA CTT CTG CTG TAA
G G G G S G G G G S S S E D D I>
-----X-----

820      830      840      850      860
      .      .      .      .      .
GAG GCC GAC CAC GTA GGC TTC TAT GGT ACA ACT GTT TAT CAG TCT CCT
CTC CCG CTG GTG CAT CCG AAG ATA CCA TGT TGA CAA ATA GTC AGA GGA
E A D H V G F Y G T T V Y Q S P>
----- I-Ad  $\alpha$ -1 domain -----

870      880      890      900      910
      .      .      .      .      .
GGA GAC ATT GGC CAG TAC ACA CAT GAA TTT GAT GGT GAT GAG TTG TTC
CCT CTG TAA CCG GTC ATG TGT GTA CTT AAA CTA CCA CTA CTC AAC AAG
G D I G Q Y T H E F D G D E L F>
-----

920      930      940      950      960
      .      .      .      .      .
TAT GTG GAC TTG GAT AAG AAG AAA ACT GTC TCG AGG CTT CCT GAG TTT
ATA CAC CTG AAC CTA TTC TTC TTT TGA CAG ACC TCC GAA GGA CTC AAA
Y V D L D K K K T V W R L P E F>
-----

970      980      990      1000      1010
      .      .      .      .      .
GGC CAA TTG ATA CTC TTT GAG CCC CAA GGT GGA CTG CAA AAC ATA GCT
CCG GTT AAC TAT GAG AAA CTC GGG GTT CCA CCT GAC GTT TTG TAT CGA
G Q L I L F E P Q G G L Q N I A>
-----

1020      1030      1040      1050
      .      .      .      .
CCA GAA AAA CAC AAC TTG GGA ATC TTG ACT AAG AGG TCA AAT TTC ACC
CGT CTT TTT GTG TTG AAC CCT TAG AAC TGA TTC TCC AGT TTA AAG TGG
A E K H N L G I L T K R S N F T>
----- I-Ad  $\alpha$ -1 domain -----

1060      1070      1080      1090      1100
      .      .      .      .      .
CCA GCT ACC AAT GAG GGT CCT CAA GCG ACT GTG TTC CCC AAG TCC CCT
GCT CGA TCG TTA CTC CGA GGA GTT CGC TGA CAC AAG GGG TTC AGG CGA
P A T N E A P Q A T V F P R S P>
-----X----- I-Ad  $\alpha$ -2 domain -----

1110      1120      1130      1140      1150
      .      .      .      .      .
GTG CTG CTG GGT CAG CCC AAC ACC CTT ATC TGC TTT GTG GAC AAC ATC
CAC GAC GAC CCA GTC GCG TTG TGG GAA TAG ACG AAA CAC CTG TTG TAG
V L L G Q P N T L I C F V D N I>
-----

```

F09020-6200000

58/58

FIGURE 29 (SHEET 4 OF 4)

```

      1160      1170      1180      1190      1200
      .         .         .         .         .
TTC CCA CCT GTG ATC AAC ATC ACA TGG CTC AGA AAT AGC AAG TCA GTC
AAG GGT GGA CAC TAG TTG TAG TGT ACC GAG TCT TTA TCG TTC AGT CAG
P  P  P  V  I  N  I  T  W  L  R  N  S  K  S  V>
-----

      1210      1220      1230      1240      1250
      .         .         .         .         .
ACA GAC GGC GTT TAT GAG ACC AGC TTC CTC GTC AAC CGT GAC CAT TCC
TGT CTG CCG CAA ATA CTC TGG TCG AAG GAG CAG TTG GCA CTG GTA ACG
T  D  G  V  Y  E  T  S  F  L  V  N  R  D  H  S>
-----

      1260      1270      1280      1290
      .         .         .         .
TTC CAC AAG CTG TCT TAT CTC ACC TTC ATC CCT TCT GAT GAT GAC ATT
AAG GTG TTC GAC AGA ATA GAG TGG AAG TAG GGA AGA CTA CTA CTG TAA
F  H  K  L  S  Y  L  T  F  I  P  S  D  D  D  I>
-----

1300      1310      1320      1330      1340
      .         .         .         .         .
TAT GAC TGC AAG GTG GAG CAC TGG GGC CTG GAG GAG CCG GTT CTG AAA
ATA CTG ACG TTC CAC CTC GTG ACC CCG GAC CTC CTC GGC CAA GAC TTT
Y  D  C  K  V  E  H  W  G  L  E  E  P  V  L  K>
----- I-Ad α-2 domain -----

      1350      1360      1370      1380
      .         .         .         .
CAC TGG GAG GAA GAA GAG TAC ATG CCG ATG GAA TGA
GTG ACC CTC CTT CTT CTC ATG TAC GGC TAC CTT ACT
H  W  E  E  E  E  Y  M  P  M  E  *>
-----X----- EE tag ----->

```

F09020-62E00060